I. OBJECTIVES

- Recognize when a reaction occurs.
- Record initial observation of solutions including pH, color, odor, viscosity, etc. Learn how to describe and document the results of mixing solutions.
- Identify 5 unknown solutions based on the results of initial observations and reactivity. Provide a brief proposal in your notebook addressing the basic experimental steps to meet these objectives. Once approved, work with your team of investigators to collect and analyze the data.

II. PROCEDURE DEVELOPMENT

In order to identify each of the 5 unknown solutions, develop a method to react each solution with each of the others for a total of 10 combinations. Document the result of each mixture, use the results to help identify each of the 5 solutions. After identification, write an equation for the reactions that occur. NOTE THAT NOT ALL COMBINATIONS WILL RESULT IN A CHEMICAL REACTION.

*For your reference, sample data tables are provided in section IV: DATA COLLECTION & CALCULATIONS

You will be provided one of the following sets of chemicals in bottles, each labeled with an unknown letter: